

**Sample Containers, Preservation
and Hold Time Requirements**

| Analyte | Container | Amt Req (mL) | Preservation ¹ | Holding Time ² |
|--|---------------------|------------------|---|---|
| <i>Microbiology</i> | | | | |
| Drinking Water Total Coliforms & Fecal Coliforms P/A | P, Sterilized | 100 | Cool If chlorinated, 0.008% Na ₂ S ₂ O ₃ Leave air space (~2.5 cm) | 30 hours |
| Non Potable Water Fecal Coliforms | P, Sterilized | 100 | Cool If chlorinated, 0.008% Na ₂ S ₂ O ₃ Leave air space (~2.5 cm) | 8 hours (6 hrs. transit with total of 8 hrs.) |
| Drinking Water HPC | P, Sterilized | 100 | Cool If chlorinated, 0.008% Na ₂ S ₂ O ₃ Leave air space (~2.5 cm) | 8 hours |
| LT2 Rule E coli Enumeration | P, Sterilized | 100 | Cool If chlorinated, 0.008% Na ₂ S ₂ O ₃ Leave air space (~2.5 cm) | 24 hrs |
| <i>Inorganic</i> | | | | |
| Acidity | P, G | 100 | Cool, No head space | 14 days |
| Alkalinity | P, G | 100 | Cool, No head space | 14 days |
| Ammonia Nitrogen | P, G | 400 | Cool, H ₂ SO ₄ to pH<2 | 28 days |
| BOD | P, G | 1000 | Cool | 48 hours |
| CBOD | P, G | 1000 | Cool | 48 hours |
| Dissolved Oxygen | G bottle and top | 300 | None | Analyze immediately |
| Chloride By IC | P, G | 50 | Cool | 28 days |
| Sulfate By IC | P, G | 50 | Cool | 28 days |
| Nitrate By IC | P, G | 50 | Cool | 48 hours |
| Nitrite By IC | P, G | 50 | Cool | 48 hours |
| Bromide By IC | P, G | 50 | Cool | 28 days |
| Fluoride By IC | P, G | 50 | Cool | 28 days |
| Chlorine, Tot.& Free Res. | P, G | 50 | None | Analyze Immediately |
| COD | P, G | 50 | Cool, H ₂ SO ₄ to pH<2 | 28 days |
| Conductivity | P, G | 100 | Cool | 28 days |
| Oil & Grease | G only | 1000 | Cool, H ₂ SO ₄ or HCl to pH<2 | 28 days |
| pH | P, G | 25 | None | Immediately analyze |
| Solids, Settleable | P, G | 1000 | Cool | 48 hours |
| Solids, TDS | Plastic | 200 | Cool | 7 days |
| Solids, TSS | Plastic | 200 | Cool | 7 days |
| Turbidity | P, G | 100 | Cool | 48 hours |
| Hex Chromium | P | 500 | Cool, (pH 9.3 – 9.7) | 24 hours |
| Hydrazine | P or G | 500 | HCl to pH < 2, Cool | 7 days |
| Total Phosphorus | P | 250 | Cool, H ₂ SO ₄ to pH<2 | 28 days |
| TOC | Amber G | 250 | Cool, HCl to pH<2 | 28 days |
| <i>Organics</i> | | | | |
| VOCs | Pre- cleaned, | 44 mLs per vial, | HCl pH<2, No head space, Cool | 14 days |

| Analyte | Container | Amt Req (mL) | Preservation ¹ | Holding Time ² |
|-------------------|--|---|---------------------------|--|
| | certified glass VOA vials pre-preserved with 25 mg Ascorbic Acid | collected in duplicate | | |
| HAA5s | Pre- cleaned, certified amber Boston round glass bottles with Teflon caps pre-preserved with 25 mg Ammonium Chloride | 250 mLs collected 40 mLs extracted | Away from light, Cool | 14 days 7 days for extract 14 days if extract preserved at – 10 Deg C or less. |
| <i>Metals</i> | | | | |
| Metals, except Hg | P, G | 100 | HNO ₃ to pH<2 | 180 days |
| Metals, Hg | P, G | 100 | HNO ₃ to pH<2 | 28 days |

¹Note: Cool = above freezing to 6°C

²Note: Hold times from sample collection

CWM Environmental, SM-001T_06 effective 01/01/15

(Replaces SM-001T_05 Effective 06/21/12)

(Replaces SM-001T_04 Effective 1/1/10)

(Replaces SM-001T_03 Effective 05/01/07, added parameters)

(Replaces SM-001T_06 Effective 11/4/15)